

# SERENU™ PROFESSIONAL

## LUMISPA iO® and LUMISPA CLEANSER

---



### MEET SERENU PROFESSIONAL LUMISPA iO AND LUMISPA CLEANSER

Experience a professional-grade cleansing ritual designed to help skin feel refreshed, balanced, and ready for the rest of your routine. LumiSpa iO pairs soft silicone treatment heads with precise Micropulse Oscillation Technology to provide a thorough yet gentle cleanse, while Serenu Professional LumiSpa Cleanser cushions the skin with a silky lather that rinses clean.

Formulated with Hibiscus Rosa-Sinensis Flower Extract, Laminaria Digitata Extract, and Alpha-Glucan Oligosaccharide, the cleanser helps comfort, hydrate, and balance skin. Together with the HydraShield Complex, these ingredients support the skin's moisture barrier and leave the complexion feeling supple and cared for.

*Serenu™ Professional is a personalized skincare system inspired by India's treasured natural ingredients to target specific skin concerns with powerful, results-driven formulas.*

---

### WHY YOU'LL LOVE IT

- Polishes skin for a smooth, soft feel.
  - Deeply cleanses while helping keep skin calm and comfortable.
  - Balances cleansing and hydration for supple-feeling skin.
  - Developed with an ingredient that helps protect the skin's microbiome.
- 

### WHAT POWERS IT

- **HydraShield Complex**—provides hydration and helps protect the moisture barrier while cleansing.
- **Hibiscus Rosa Sinensis Flower Extract and Laminaria Digitata Extracts**—deliver skin conditioning benefits to help keep the skin soft and supple.
- **Alpha-Glucan Oligosaccharide**—a prebiotic that contributes to balanced, comfortable skin.

## HOW TO USE IT

Moisten face with warm water and apply an ample amount of Serenu Professional LumiSpa Cleanser evenly over the face, avoiding the eye area. If contact occurs, rinse thoroughly with water. Use with Serenu Professional LumiSpa® iO as directed in the device manual.

## TRY THESE WITH IT

Serenu™ Professional Brightening Toner  
Serenu™ Professional Galvanic Spa®  
Serenu™ Professional Brightening Serum  
Serenu™ Professional Gel Cream

---

## LEARN MORE ABOUT IT

### **What makes Serenu Professional LumiSpa iO different from other cleansing devices?**

Serenu Professional LumiSpa iO features proprietary Micropulse Oscillation Technology combined with a soft silicone head, enabling a thorough yet comfortable cleanse. Unlike traditional brushes, it is microbiome-friendly, respects skin's pH, and helps preserve its natural moisture barrier while providing a refreshed, spa-like experience.

### **Do I need to use Serenu Professional LumiSpa Cleanser with Serenu Professional LumiSpa iO, or can I use other cleansers?**

Serenu Professional LumiSpa Cleanser was formulated to work with the Serenu Professional LumiSpa iO device to deliver the best results.

### **What makes Serenu Professional LumiSpa Cleanser and Serenu Professional LumiSpa iO unique?**

Together, they form the first step in the Serenu Professional system. The Serenu Professional LumiSpa iO device provides a thorough yet gentle cleanse, while Serenu Professional LumiSpa Cleanser cushions and conditions the skin, removing impurities and leaving it feeling refreshed and hydrated.

### **What makes the Serenu Professional LumiSpa iO Treatment Head unique?**

The soft, nonabrasive silicone surface contains antimicrobial silver and is more hygienic than brushes commonly used with skincare devices. The sophisticated design helps to maximize product efficacy while making skin more receptive to subsequent products.

### **How long will my Serenu Professional LumiSpa iO Treatment Head last?**

Each treatment head is designed to last for approximately 180 professional treatment sessions. Replacement frequency may vary depending on usage patterns in a salon environment.

**Who should use the Serenu Professional system?**

This system is designed for anyone who wants a professional-inspired skincare experience at home. It is suitable for daily use and helps maintain balanced, healthy-looking skin. Whether you're seeking a comfortable daily cleanse or a complete routine, the Serenu Professional system offers a holistic approach.

**Which aspects of Serenu Professional LumiSpa Cleanser are drawn from India's treasured natural ingredients?**

The Serenu Professional LumiSpa Cleanser incorporates traditional botanicals such as Hibiscus Rosa-Sinensis Flower Extract and Laminaria Digitata Extract. These ingredients are traditionally valued for their conditioning and balancing properties, aligning with the heritage of holistic beauty.

**How long does Serenu Professional LumiSpa iO last after a full charge with the magnetic charger?**

Serenu Professional LumiSpa iO typically allows up to 50 sessions per full charge. Usage time may vary depending on frequency and duration of treatments.

**What makes Serenu Professional LumiSpa iO different from other cleansing devices?**

Serenu Professional LumiSpa iO features proprietary Micropulse Oscillation Technology combined with a soft silicone head and the specially formulated Serenu Professional LumiSpa Cleanser. Together, they provide a thorough yet comfortable cleanse. Unlike traditional brushes, the system is microbiome-friendly, respects skin's pH, and helps preserve its natural moisture barrier while leaving skin feeling refreshed and balanced.

**Do I need to wash my makeup off before using Serenu Professional LumiSpa iO?**

Serenu Professional LumiSpa iO conveniently and effectively cleanses away makeup, so you do not need to wash makeup off before using the device. However, some people may like to double cleanse to remove heavy makeup. In that case, you may want to use another cleanser or use your Serenu Professional LumiSpa iO twice. Please note that Serenu Professional LumiSpa Cleanser should not be used on the eye area, so you should use a separate eye makeup remover to remove makeup on and around the eyes.

**What other products are in the Serenu Professional system?**

The complete system includes Serenu Professional Brightening Toner, Gel Cream, Galvanic Spa and Brightening Serum—each designed to work with Serenu Professional LumiSpa iO and Cleanser for professional-inspired results.

**How do I know if I am using the right motion and right amount of pressure?**

It's best to move the device across your skin in slow, broad strokes. The device can prompt you in using the proper motion and right amount of pressure. If you use a scrubbing motion, the treatment head will pause and vibrate three times before resuming. If you press too hard, the treatment will pause, followed by one long vibration before resuming. Depending on your routine, you will receive a prompt to move from one area of the face to the next.

**Where can I use my Serenu Professional LumiSpa iO?**

Serenu Professional LumiSpa iO is designed for use on the face, neck, and décolleté. We recommend using it only on these areas as part of your daily skincare routine.

**Can I use Serenu Professional LumiSpa iO in the shower?**

Yes. Serenu Professional LumiSpa iO is waterproof and safe to use in the shower. To maintain your device, avoid dropping it or leaving it submerged in water.

---

**PRECAUTIONS**

For cosmetic use only. Avoid direct contact with eyes. Discontinue use if irritation occurs. Individuals with sensitivities should carefully review ingredient listings prior to use.

---

**WHAT'S IN IT**

Aqua (Water), Glycerin, Caprylyl/Capryl Glucoside, Ammonium Polyacryloyldimethyl Taurate, Glyceryl Stearate, Caprylic/Capric Triglyceride, Isoamyl Laurate, Panthenol, Potassium Cetyl Phosphate, Polyacrylate Crosspolymer-6, Hibiscus Rosa-Sinensis Flower Extract, Alpha-Glucan Oligosaccharide, Ectoin, Disodium EDTA, T-Butyl Alcohol, Sodium Hyaluronate, Sodium Carboxymethyl Beta-Glucan, Laminaria Digitata Extract, 1,2-Hexanediol, Sodium Hydroxide, Glycolic Acid, Hydrogenated Lecithin, Ceramide 1, Ceramide 2, Ceramide 3, Ceramide 4, Ceramide 6 II, Ethylhexylglycerin, Cetyl-PG Hydroxyethyl Palmitamide, Phenoxyethanol, Cholesterol, Chloroacetic Acid, Sodium Benzoate, Nitrogen.